


Documentation to HPLC Data Vision

Institute	University of Chemistry and technology, Prague	
Faculty	Faculty of Chemical Engineering	
Department	Department of Mathematics, Informatics and Cybernetics	
Author	Jakub Vencel	
Contact	jakub.vencel24@gmail.com	
Date	May 1, 2025	
App Version	1.0	

1 About

This application provides users with the ability to interact with data from High Performance Liquid Chromatography. Its features include extracting data from *.txt files, exporting data to Excel sheets, generating L^AT_EX code for PDF protocols (such as the template described in section 3), and displaying chromatograms in a separate window.

1.1 Technical Information

1.1.1 Compatibility

Application is compatible with LabSolutions (ver. 5.99). It is specifically designed for HPLC with PDA detector.

1.1.2 System Requirements

Operating system: Windows

Hardware: unspecified


1.1.3 Source Codes & Architecture Diagram

Programming

Program was written in Python 3.13. Used modules:

1. CustomTkinter [1]
2. TeX Live [6]
3. OpenPyXL [4]
4. NumPy [3]
5. Matplotlib [2]
6. Pillow [5]

Figure 1: Data flow diagram



1.1.4 Support

For further information contact me through my e-mail: jakub.vencel24@gmail.com.

2 How to use

2.1 Instructions

1. Export files from LabSolutions (ver. 5.99) to *.txt; use default item separator (tabulator).
2. Load HPLC Data Vision (Figure 2) and select folder, where you store exported *.txt files. (Warning: Selected folder should contain just *.txt from LabSolutions (ver. 5.99)!)
3. Mark off desired operations (Figure 3).
4. Loading...
5. In a new window with chromatogram (Figure 4) you can choose wavelength and displayed file.

Figure 2: Main Window




Figure 3: Select Parameters





Figure 4: Chromatogram Viewer



2.2 Disclaimer

Not following instructions mentioned above may lead to unpredictable behavior!

3 Protocol Template


3.1 Measurement Information

File Name
Sample Name
Date
Time
Measuarment time (min)
Method File
Data File
Batch File
Instrument
Detector

3.2 Peaks

peak id	τ min	area	height	c
------------	---------------	------	--------	-----

3.3 Chromatogram



Contents

1	About	1
1.1	Technical Information	1
1.1.1	Compatibility	1
1.1.2	System Requirements	1
1.1.3	Source Codes & Architecture Diagram	1
1.1.4	Support	2
2	How to use	2
2.1	Instructions	2
2.2	Disclaimer	3
3	Protocol Template	4
3.1	Measurement Information	4
3.2	Peaks	4
3.3	Chromatogram	4
A	Acknowledgment	6

A Acknowledgment

References

- [1] Customtkinter. <https://customtkinter.tomschimansky.com/>.
- [2] Matplotlib. <https://matplotlib.org/>.
- [3] Numpy. <https://numpy.org/>.
- [4] Openpyxl. <https://openpyxl.readthedocs.io/en/stable/>.
- [5] Pillow. <https://pillow.readthedocs.io/en/stable/>.
- [6] Tex live. <https://www.tug.org/texlive/>.